

enesers

Operating Lamps Sinar 7050



Specifications

Model	Sinar 70	Sinar 50
Luminance(Lux)	160000	120000
Diameter of Light Head(CM)	70	50
Number of Bulbs	108	54
Color Temperature (K)	4300±200	4300±200
Color Render Index (Ra)	95	95
Red Color Index (R9)	94	94
Light Field Diameter (CM)	12~35	12~35
Illumination Depth (CM)	>144	>145
Working Distance (CM)	50~180	50~180
Bulb Service Life (H)	>50000	>50000
Weight	75 Kg	
Loading Capacity of Ceiling	More than 500 Kgs	
Kamera	1/3 CMOS Sensor	
Monitor	22 inch	



KEMENKES RI
AKD 21603910692

Manufactured By:

enesers

Taman Tekno BSD City Sektor XI

Blok K-2 No.12

Setu, Tangerang Selatan - Banten 15314

Telp : +6221 7566360

Fax : +6221 7566361

Email : info@enesers.com

Website: www.enesers.com

Distributed By:

PT.PRIMA CITRA INOVINDO Cibubur
Times Square, Blok B1 No,17 Jl.Alternatif
Cibubur KM 3. Kota Bekasi Jawa
Barat.17435Tlep : 021-50207596
Email : info@pci-tech.co.id
Web : pci-tech.co.id